

Docker Course

Bioinformatics software tests in docker containers

- Alignment and quantification tools for RNA sequencing data -

DTU Biosustain
Data Science Platform


Sebastian Schulz

Senior Data Scientist

sebschu@dtu.dk

26 November 2025

Major bioinformatic steps in RNA sequencing (RNAseq)




Major steps in RNAseq processing

- Alignment: map reads to a reference sequence (e.g. genome)
- Read quantification: count how many reads are aligned to a specific genomic region (feature)

One major downstream analysis goal

- Differential expression (DE) analysis comparing read counts of a specific genomic feature between biological experiments


Healthy cells



Cancer cells




Data processing – nf-core/rnaseq pipeline



Pro- and eukaryotic data processing with the nf-core/rnaseq pipeline

- The pipeline **succeeds with eukaryotic data sets** (e.g. mouse)
- The pipeline **fails with prokaryotic data sets** (e.g. *E. coli*) due to
 - Suboptimal quality of bacterial gtf annotation files and therefore an incompatibility issue of STAR and Salmon
 - Successfully addressed by modifying the gtf file and pipeline parametrization ([github issue #1512](#))
- The aligner-quantifier tool set, STAR + Salmon, is **not necessarily the first choice for the processing of prokaryotic data set**
 - STAR was developed for read alignment using eukaryotic data sets
- Is the combination of STAR + Salmon still applicable to prokaryotic data sets?
 - How do the results obtained from STAR + Salmon compare to those from tool sets that are (commonly) used in prokaryotic transcriptomics, like Bowtie2 + featureCounts?




<https://nf-co.re/rnaseq/3.16.1>

Alignment and quantification software

Read alignment

	Splice-junction detection	Application
STAR Dobin <i>et al.</i> , 2012	yes	<i>De novo</i> splice junction detection (eukaryotes)
Bowtie/ Bowtie2 Langmead <i>et al.</i> , 2009 Langmead and Salzberg, 2012	no	Established aligner for prokaryotes (also for eukaryotes)



MMP: Maximal Mappable Prefix
From Dobin *et al.*, 2012

Read quantification

	Quantification model	Application
Salmon Patro <i>et al.</i> , 2017	Complex statistical model considering 5' and 3' bias, GC-content bias, etc	More sophisticated quantification model for more precise results considering different biases
featureCounts Liao <i>et al.</i> , 2014	Conceptually simpler counting mechanism	Performs well on one-isoform genes (Parelo <i>et al.</i> , 2014; tested on eukaryotic data sets), which would be desirable for prokaryotic data

Dobin *et al.*, 2012, <https://doi.org/10.1093/bioinformatics/bts635> (STAR)
 Langmead *et al.*, 2009, <https://doi.org/10.1186/gb-2009-10-3-r25> (Bowtie)
 Langmead and Salzberg, 2012, <https://doi.org/10.1038/nmeth.1923> (Bowtie2)
 Liao *et al.*, 2014, <https://doi.org/10.1093/bioinformatics/btt656> (featureCounts)
 Patro *et al.*, 2017, <https://doi.org/10.1038/nmeth.4197> (Salmon)
 Parelo *et al.*, 2024, <https://doi.org/10.1093/nargab/lqae020>

Comparison of alignment and quantification software

Goal

- To test how different aligner-quantifier combinations perform on prokaryotic data sets **outside of the nf-core/rnaseq pipeline**
 - Example: Does Bowtie2 + featureCounts perform similar or different compared to STAR + Salmon?

Approach

- Substituting software tools directly in the nf-core/rnaseq pipeline just for testing purposes would be too labour-intensive and time-consuming
- Test different aligner-quantifier combinations **in Docker containers** to ensure **maximal reproducibility** for us and most importantly when **shared with the nf-core community**
 - As these analysis should form the basis for the potential further development of the nf-core/rnaseq pipeline, we aimed at ensuring a high-level of standardization in our analyses

Quantifier		
Aligner	Salmon	featureCounts
	STAR	Pipeline default
	Bowtie	
	Bowtie2	This Docker course

Reduced version of entire test plan




ACKNOWLEDGEMENTS

**novonordisk
fonden**

Supplementary slides


Create the codespace

- Log into your github account
- Open the training repo: https://github.com/biosustain/dsp_transcriptomics_training
- Create codespace **with 4 cores** which we will use to run code




Select **4-cores** for machine type before creating the codespace

Create codespace for
biosustain/dsp_transcriptomics_training




1. Data processing - Read alignment with STAR

- Spliced Transcripts Alignment to a Reference (STAR)
- Very fast aligner
- Alignment against genome (or transcriptome)
- *De novo* detection of splice junctions (no prior annotation of splicing event required)



Detection of splice junctions with STAR



MMP: Maximal Mappable Prefix

1. Data processing - Read quantification with Salmon

- More complex quantification procedure compared to other read counters/quantifiers
- Handles counting of multi-mapping reads
- More accurate quantification of reads achieved by considering sample-specific parameters and biases of RNAseq data:
 - positional biases in coverage
 - sequence-specific biases at the 5' and 3' ends of sequenced fragments
 - fragment-level GC bias
 - strand-specific protocols
 - fragment length distribution

